

# Product Information

## Digital Blood Pressure Monitor



**Brand :** Health Pro

**Model :** HP-ES004

**Origin:** PRC

**Feature:**

- Australian Brand
- 99\*2 sets memory for 2 users
- English voice
- Date & time setting
- Delete memory function
- Oscillometric measuring method
- Screen size: 55.2\*75\*2.8mm

**Size :** 136mm(L) x 102mm(W) x 73mm(H)

**N.W.:** ~236g (without batteries)

**Certificate :** CE

**Specs:**

- Operating Voltage: DC6V (AA Batteries)
- Display: LCD Digital Display
- Inflatable method: Automatic Inflation
- Resolving power: 1mmHg/0.1kPa
- Cuff size: 22 ~ 32 cm
- Operating Temperature: 5-40°C
- Operating Humidity: ≤5%-80%RH
- Storage Temperature: -20 - +55°C
- Storage Humidity: ≤93%RH
- Pressure range: 0-37.3kpa(0-280mmhg)
- Pulse range: 40-199 beats/min
- Accessories Included: Wrist Band, Product Manual, Storage Bag