

Bedhead Trunking

The Bedhead Trunking is specially designed to fit Terminal Units of British standard and other standards.




To comply with the international safety standards to separate electricity supply to gas supply, the Bedhead Trunking consists of individual ducts to accommodate high voltage wiring, low voltage wiring and a maximum of 4 different gas piping.


At customer's need, an internal reading light can be fitted in the bottom of the Bedhead Trunking. Facial plates are mounted by interlocking joint enable easy removal for servicing.

Specification: 1.2m length aluminium alloy bedhead trunking, with gas outlet (O2, VAC, AIR) x3, nurse call x1, emergency call x2, electric socket x4, Room light & Reading light with with ON/OFF switch.

SBDA-B



SBDA-K



SBDA-R

