


Central Control Panel 6 Compartment



Function:

- Clock
- Timer
- Gas alarm
- Intercom / Telephone
- Temperature control / Humidity display
- Light control

Model: HP-OR-TCP

These systems constitute the information monitoring of the operating room. Every part could be displayed independently by the menu option buttons, in addition to the clock timing system, the rest of the monitoring system are connected by a computer via RS232 port to the underlying hardware interface, by the hardware with computer real-time and high-speed communication, data exchanged to the underlying input/output control electronic components.