



Product Information

Fingertip Pulse Oximeter (HP-ES005)



Technical Specification:

- CE Certificate
- Automatic shutdown
- Fingertip: Recommend left middle finger; fingernail point upward
- Colored TFT LCD display
- Measurement Range: From 0% to 100%
- Oxygen measurement accuracy in the range of 70% - 100%: $\pm 2\%$
- Heart rate measurement range: 30 to 250 BPM
- Heart rate accuracy: ± 2 BPM
- Weight: 50g (Included battery)
- Power supply: 2 x AAA (1.5V Alkaline battery)
- Display dimension: 22mm x 11mm
- Device dimension: 57 x 31 x 32mm

Goodwin Health Care Ltd. 嘉榮醫療有限公司

Tel: (852) 2408 8838

Fax: (852) 2408 1329

For WhatsApp Only: (852) 9872 0697

Email: hksales@goodwin.com.hk

Website: www.goodwin.com.hk

