

Surgical LED Light

Operating LED Light, medical equipment double head lamp CreLed 3300/3300: Double Dome Operating Room Light, the illumination can be adjusted freely in the range of 25%-100%, which can relieve the fatigue of the eyes of doctors who have been operated on for a long time. Double Dome Surgery Room Light, the digital memory function automatically records the proper station manes and does not need to be adjusted again. Double Dome Operating Theatre Light, in the condition of clear and comfortable, you can see the details of each exposure group and the details in the shadow of each tiny red shadow.



Model No.: CreLed 3300/3300

Type: Operating Light	Color Temperature(K):3500-5000
Kind: Ophthalmology	Lighting Depth(mm):≥1200mm
LCD Display: With LCD Display	Electric Spot Size Adjustment: Yes
Source: Air	Mode: Background, Surface, Deep Cavity, Normal
Anesthesia: Semi Open Suction	Color Rendering Index(Ra): 95
Type of Operating Light: Double Operating Light	Diameter Of Light Head: 600 mm/ 600mm
LED Life: More Than 50000 Hrs	Quantity Of LED Bulbs: 45/45
Spot Diameter(mm):150-280	Illumination: 60000-130000 LUX