



Health Pro

Unique Innovation with
All Rights Reserved

Developing of ward and patient care



1st Generation: Movable Gas Cylinder



Old fashioned gas cylinder is unsafe. Time-consuming and troublesome.

2nd Generation: Bed Head Unit



- One straight line design bring patients uneasy feelings.
- Exposed gas outlets and power sockets may cause safety risks when do cleaning.
- Aluminum alloy housing is easy to conduct electrostatic or cause current leakage.
- Unable support to digital information which prevent modern clinical development and upgrading.

3rd Generation: Integrated Ward Care System



- To meet demands of modern clinical development, Health Pro Pty Ltd developed Integrated Ward Care System which bring hospital's management to a new stage.
- The system extends digital information to each bed. It extremely improves the efficiency and relieves medical staff from the most tedious work of patient caring.
 - Graceful soft gentle line design help patient release stress feelings.
 - Well partitioned and cover protected gas outlet and power socket for more safety concern.
 - The smooth and bold cabinet design offer enough storage space for patient belongs.

Perfect Integration of Different Functions
Caring Service for Patients and Medical Staff



Integrated Panel Series

Integration of medical gas outlet, power socket, nurse call, internet and etc.



Body made by high strength engineering plastic in one piece. Compact in structure, simple and smooth in design

Equip with leakage protection switch, 4pcs power sockets and 1pc grounding terminal for safety connection of the medical equipments

concealed built-in high and low current partitioned by different segment and channel to ensure safety

Reserved place for internet and other low current

Gas outlet with movable shield for dustproof and easy maintenance



Normal Integrated Panel



Support to traditional bed head information card

Smart Integrated Panel



Smart digital information system and nurse call system supported by PAD on the panel which compatible to mobile medical system and smart home medical system interface, support to remote consultation.

Personal Stuff Storage Cabinet

The similar concept of aircraft luggage compartment near the patient for easy access. Simple and bold design offer enough storage space for patient belongs. It not only improves the effective utilization of the ward space, but also solves the difficulty of the storage management. It provides a feasible solution for the hospital ward environment optimization.

Body made by high strength engineering plastic in one piece in compact structure. Simple, smooth and grace in design upgrade ward environment. Multi-function integrated and easy to control.

Electronic lock function provides security for patient's personal belongs.

Built-in UV light for disinfection



A strip LED light is designed at the bottom of the cabinet

Integrated Ward Care System Smart 2beds



Each bed equip with 2 storage cabinet, 1 set smart Integrated pannel and digital nurse call system

Integrated Ward Care System-Normal 2beds



Each bed equip with 2 storage cabinet, 1 set normal inegrated pannel

Integrated Ward Care System-Smart VIP



Integrated Ward Care System-Normal VIP



Each bed equip with 3 storage cabinet, 1 set smart integrated panel and digital nurse call system

Each bed equip with 3 storage cabinet, 1 set normal integrated panel

Ward Area Smart Digital Management and Nurse Call System Diagram

