




LED Backlight X-Ray Viewer Box (2 Bay)

Features

- ◆ **LED Light**, bright immediately, bright and with 20000 hours life
- ◆ **Special Diffuser**, making a regular luminous density, and the accordant luminance from the 4 ends to the center of the X-ray Viewer Box house
- ◆ **Colour Temperatures:6200K**
- ◆ **Glare-free Reading**, 4 shutters around the X-ray Viewer Box house are used to adjust the Image Area, which can be diminished from the full Image Area to the minimum 60 x 60mm., but with shape-edge collimation even of small Image Areas, so that the doctor can read more carefully and make the right judgment
- ◆ **Durable House**, five layers of fancy electrostatic spraying, making the house anti-scratch, abrasion resistant and anti-corrosion, furthermore, the outer layer of the house is made of high strength tempered glass while the inter layer is made of high quality rubber sheet
- ◆ **Types:** Wall-mounted or Desktop
- ◆ **Bay:** 2 bays

Model: HP-OR-XB2

Size: 1000(w) x 625(h) x 102(d)





LED Backlight X-Ray Viewer Box (3 Bay)

Features

- ◆ **LED Light**, bright immediately, bright and with 20000 hours life
- ◆ **Special Diffuser**, making a regular luminous density, and the accordant luminance from the 4 ends to the center of the X-ray Viewer Box house
- ◆ **Colour Temperatures:6200K**
- ◆ **Glare-free Reading**, 4 shutters around the X-ray Viewer Box house are used to adjust the Image Area, which can be diminished from the full Image Area to the minimum 60 x 60mm., but with shape-edge collimation even of small Image Areas, so that the doctor can read more carefully and make the right judgment
- ◆ **Durable House**, five layers of fancy electrostatic spraying, making the house anti-scratch, abrasion resistant and anti-corrosion, furthermore, the outer layer of the house is made of high strength tempered glass while the inter layer is made of high quality rubber sheet
- ◆ **Types:** Wall-mounted or Desktop
- ◆ **Bay:** 3 bays

Model: HP-OR-XB3

Size: 1300(w) x 630(h) x 102(d) mm

