



Mobile Negative Pressure Isolation Ward

This is our new design of isolation ward with a short video share to you. Our R & D active develop this type. We think it suitable to many countries under the epidemic situation. It's size same as a 20 feet container. Easy for transportation by truck.

Mobile Negative Pressure Isolation Ward overall size 6m(L)*2.4m(W)*2.8m(H) same as a 20 feet container easy for transportation by container trunk remove to other virus outbreak city when they haven't enough Isolation Ward.

No need extra AHU. This is integration design one set included Air condition, pressurized fan, exhaust fan, Silver Ion H14 grade HEPA filter, air compressor, vacuum pump of medical gases for bedhead Trunking, on site only need have a Oxygen cylinder for supply oxygen. Interlock Pass through window, automatic doors.

On site after fixing just connect the power, water supply and drain can running for use immediately.

Mobile negative pressure Isolation ward use for once Virus outbreak district also the neighboring hospital haven't or not enough Isolation ward can fastest take action send these wards with very short time to the place saving people also effectively make the infection control the virus.



Download link of Health Pro Model V1 - Mobile Negative Pressure Isolation Ward remote video.

https://drive.google.com/file/d/1CWQzBKyxJ7aWaEpGR3ZVjg8Sb3_VizY9/view?usp=sharing

